

SCHURE SPINE

#800-0235



Adjust
intraoperatively via
the crank handle to
desired lordosis to get
optimal surgical site
exposure.



Adjusts
intraoperatively to
the desired lordosis
for optimal surgical
site exposure




A secure, comfortable and adjustable spinal platform

The Schure™ Spine Frame provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions, and disc surgery. Tapered adjustable pads provide optimal patient safety and comfort.

The Schure Spine Frame has improved features over the current Wilson™ Frame, has a patient weight capacity of 500 lbs. (227 kg) and is C-arm compatible. It can be used with any O.R. table, and fits the Jackson Spine Table and Allen Flex Frame as well.

Features

- ◆ Adjusts intra-operatively to the desired lordosis to get optimal surgical site exposure
- ◆ 500 lb patient weight capacity
- ◆ Frame comes complete with Schure™ Spine Frame Dolly, Crank Handle, Disposable Headrest, Spine Frame Pad Set, and 2-Arm Cradles
- ◆ Radiolucent and allows for unrestricted C-arm access
- ◆ Can be used with any O.R. table and fits the Jackson Spine Table and Allen Flex Frame



#800-0235 Schure Spine Frame

Replacement Products

#508-0472 Schure Spine Frame Pad Set
#800-0277 Schure Spine Frame Crank Handle

Optional Accessories

#508-0373 Disposable Laminectomy Arm Cradles, 12 prs/cs
#508-1340 Disposable Richard Slotted Head Rests, Extended, 6/cs
#800-0074-SF Schure Spine Frame Dolly